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
Figure 1. T- cells isolated from the spleen were fluorescently stained with DiR and i. v. injected (5×10^6 cells/ mouse) into a Nu/ Nu mouse. Images above taken 24hrs post injection with IVIS Spectrum show cells homing to the spleen.

In Vivo Imaging of DiR Stained Spleen T-cell Distribution


Figure 2 . DiR stock was prepared by dissolving 25 mg in 3 mL ethanol . Working solution of 320 pg/ mL was prepared by diluting 199 μ L of stock solution in 5 mL PBS. T- cells isolated from the spleen were incubated with 320 pg/ mL DiR. After 30 min incubation, cells were spun down for 3 min at 1000 rpm at 4 $^{\circ}$ C resulting in a blue pellet . Cells were washed twice in PBS and injected intravenously (5×10^6 cells/ mouse). Control group was injected with 5×10^6 cells/ mouse in PBS. Mice were imaged with IVIS Spectrum at 10 min, 1hr, 6hr and 24 hrs post injection . Ideal filter set for DiR imaging is 710 nm excitation and 760 nm emission . Mice were imaged dorsally as well as ventrally at all time points . Brain, bones, spleen, liver, lungs and kidneys were harvested for ex vivo imaging 24 hrs post injection .

Non- invasive in vivo imaging showed the homing process of injected T cells to the liver and spleen in real time, which was confirmed by ex vivo imaging .